



Advanced X-ray Characterization of Hierarchical Systems: from Biotissues to Electronic Devices.

Davide Altamura,^{a,*} Francesco Scattarella,^a Alberta Terzi,^a Teresa Sibillano,^a Rocco Lassandro,^a Liberato De Caro,^a Dritan Siliqi,^a Massimo Ladisa,^a Cinzia Giannini.^a

Label-free Multiscale X-ray Imaging can be performed by Table-Top equipment, exploiting both Absorption (usually referred to as X-ray Transmission Microscopy “XTM”) and Scattering contrast. X-ray scattering techniques can provide a large amount of structural and morphological information, both at the atomic and nano-scale, and are thus particularly suited to study composite/nanostructured materials. The crystalline components are mainly studied by Wide Angle X-ray Scattering (WAXS), providing information on the crystallinity and crystalline phases, as well as on possible texture. The nanoscale structure/morphology can be assessed based on the Small Angle X-ray Scattering (SAXS) signal, and related to the possible crystallinity through combined SAXS/WAXS mapping. SAXS/WAXS/XT Microscopy is particularly suited to study biological tissues (or any structured material) with nano and/or atomic scale periodicity, such as calcified healthy or pathological tissues.¹ In the biomedical field, this approach allows for example to study mineralized bio-scaffolds, or reveal osteogenic differentiation of stem cells through the analysis of nanocrystalline differentiation products.² The availability of high brilliance X-ray micro-sources for laboratory equipment allows nowadays to perform the aforesaid advanced X-ray characterization, in both transmission and reflection geometries, in the home laboratory, being such X-ray sources considered as “synchrotron-class”. Moreover, suitable data treatment can further enhance the performances of the experimental equipment, returning in several cases results comparable to those obtained at a synchrotron beamline. A in-house developed package (SUNBIM) for the collection and analysis of X-ray microscopies with absorption and/or diffraction contrast, as well as data reduction for transmission and reflection geometries is freely available at <http://www.ba.ic.cnr.it/softwareic/sunbimweb/>, and constantly updated. The synchrotron-class micro-source combined with the SAXS/WAXS system installed at IC-Bari (XMI-L@b) has been successfully applied to the study of free-standing bio- and composite materials, as well as to polymer nanofibers for optoelectronic devices. Moreover, it has been successfully applied in grazing incidence reflection geometry (GISAXS/-GIWAXS) for the study of several nanostructured films for advanced electronic applications, in particular based on 2D nanocrystals.³

References

- 1 C. Giannini, M. Ladisa, V. Lutz-Bueno, A. Terzi, M. Ramella, L. Fusaro, D. Altamura, D. Siliqi, T. Sibillano, A. Diaz, F. Boccafoschi, O. Bunk, X-ray scanning microscopies of microcalcifications in abdominal aortic and popliteal artery aneurysms, *IUCrJ* 6 (2) (2019) 267–276. doi:10.1107/S2052252519001544.
- 2 A. Procopio, E. Malucelli, P. Alexandra, C. Cappadone, G. Farruggia, A. Sargenti, S. Castiglioni, D. Altamura, A. Sorrentino, C. Giannini, E. Pereiro, P. Cloetens, J. A. M. Maier, S. Iotti, Chemical Fingerprint of Zn–Hydroxyapatite in the Early Stages of Osteogenic Differentiation, *ACS Cent. Sci.* 5 (8) (2019) 1449–1460. doi:10.1021/acscentsci.9b00509.
- 3 F. Scattarella, E. Altamura, P. Albanese, D. Siliqi, M. Ladisa, F. Mavelli, C. Giannini, D. Altamura, Table-top combined scanning x-ray small angle scattering and transmission microscopies of lipid vesicles dispersed in free-standing gel, *RSC Advances* 11 (1) (2021) 484–492. doi:10.1039/D0RA08581B.

^a CNR - Istituto di Cristallografia, Via Amendola 122/O, 70126 Bari, Italy

Creative Commons Attribuzione - Non commerciale - Condividi allo stesso modo 4.0 Internazionale

† oral communication at 1 st Conference on Crystallography, Structural Chemistry and Biosystems, (Catania) 04-06/10/2021